


Eakin Ring

- Indications
 - Seal or barrier for uneven skin around ostomies, wounds, and/or fistulas
- Open package and remove release paper
- Pinch/flatten outer edges of ring to create a lower profile




Eakin Ring

- Stretch the ring to fit against the stoma
 - The Eakin Ring can/should touch the stoma
- As needed, mold/compress ring to fit down into any uneven surface of the skin



Eakin Ring

- Smooth the ring on the skin and toward the stoma



Eakin Ring

- For small, pediatric stomas, use a strip of Eakin
- Roll/mold strip of Eakin into narrow “worm”
- Mold “worm” of Eakin to fit around stoma and follow previous stepladder

